

Arthroscopy Single use Disposable Shaver Blades and Burrs

High Precision Performance Unmatched Versatility Superior Design



Linvatec

Disposable Arthroscopy Shaver Blades and Burrs

Shark Cut Linvatec Type





	TECHCORD#	Radius	Linvatec#
	T2730	3.0 mm	
Full Radius	T2735	3.5 mm	C9248
	T2742	4.2 mm	C9247A
Aggressive	T1735	3,5 mm	C9253
Aggressive	T1742	4.2 mm	C9254
Bevel Edge Cutter	T3730	3,0 mm	
Dissector	T3735	3.5 mm	C9240
Dissector	T3742	4,2 mm	C9242
Oval Burr	T1750	4.0 mm	H9103RH
Round Burr	T2743	4.0 mm	н9111, н9118

Arthrex

	The state of the s	TECHCORD#	Radius	Arthrex#
Full Radius	P	T2842	4.0 mm	AR-8400BC
Aggressive		T1842	4.0 mm	AR-8400EX
Oval Burr	P	T1850	4.0 mm	AR-74000BE
Round Burr	Co.	T2843	4.0 mm	AR-7400RBE

Stryker

Disposable Arthroscopy Shaver Blades and Burrs

Shark Cut Stryker Type





			Shark Cuty		
		TECHCORD#	Radius	Stryker#	
		T2635	3.5 mm	375-532-000	
Full Radius	Full Radius	T2642	4.2 mm	375-542-000	
Aggregative		T1635	3.5 mm	375-535-000	
Aggressive		T1642	4.2 mm	375-545-000	
Payal Edga Cuttor		T3635	3.5 mm	375-534-000	
Bevel Edge Cutter Dissector		T3642	4.2 mm	375-544-000	
Oval Burr	P	T1650	4.0 mm	375-951-112	

Round Burr



T2643

4.0 mm

375-940-012



Smith & Nephew

Disposable Artroscopy Shaver Blades and Burrs

Shark Cut Smith & Nephew Type





	TECHCORD#	Radius	Smith&Nephew#
·	T2930	3.0 mm	
Full Radius	T2935	3.5 mm	7205305
	T2942	4,2 mm	7205306
Aggressive	T1935	3.5 mm	7205345
	T1942	4.2 mm	7205459
	T3930	3,0 mm	
Bevel Edge Cutter Dissector	Т3935	3.5 mm	
	T3942	4.2 mm	7205310
Oval Burr	T1950	4.0 mm	7205326
		1	
Round Burr	T2943	4.0 mm	7205324

HRC 52 Degree Enhance the hardness of Blade tip

- Patented manufacturing method for SharkCut Blades
- Reduced Clogging
- Pure cut with less tissue tearing



TECHCORD



COMPETITER1



COMPETITER2

Techcord offers a wide range of high quality, single use Disposable Shaver Blades and Burrs for Stryker, Linvatec, Arthrex, Smith&Nephew arthroscopic surgical procedures.

Made with superiror materials, manufacturing Processes and workmanship, these products offer uncompromising performance and efficiency during tissue resection and Burring applications.

The Disposable Shaver Blades and Burrs can be used in both the large and small Shaver hand pieces and can be inserted with the cutting window in the up or down position.

Threads molded into the Blade and Burr hubs allow them to lock onto bridge cannulas for increased maneuverability and fluid management.





Disposable Small Joint and Hip Shaver Blades and Burrs

Made with the same attention to detail as their standard counterparts, the disposable small joint and hip Shaver Blades and Burrs provide outstanding performance in bone and tissue resection applications.

The disposable small joint and hip Shaver Blades and Burrs can be used in standard sized foot control and hand control Shaver hand pieces and can be inserted with the cutting window in the up or down position.



Precision It's Your Ch

(주)**텍코드** Techcord Co., Ltd.

(34014)대전광역시 유성구 테크노 1로 62-8 Tel. (042)221-0015 Fax. (042)221-0016

62-8 Techno 1-ro, Yuseong-gu, Daejeon, 34014, Republic of Korea

Tel, +82-42-221-0015 Fax, +82-42-221-0016 www.techcord.co.kr

Distributed